

ROHO Padded Seat



- The ROHO Padded Seat protects skin while sitting on a commode chair. The cushion comes standard to fit commode seats 41 cm x 41 cm (16 x 16, 9 x 9 cells) and 45 cm x 45 cm (18 x 18, 10 x 10 cells)
- Designed to fit shower chairs, the Padded Seat cushion offers the therapeutic benefits of the LOW PROFILE® cushion
- Made from high quality neoprene, the ROHO padded seat comes in a number of standard sizes, or can be custom made to suit different seat styles
- Approximate Weight: 0.8 kg
- Custom sizes and shapes are available on request

Warranty Information

- 1 Year

permobil

ROHO



Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_093330084_0103_1_2

Mobile Shower Commode –

Composite: 05_413_0103_1_2

Click [here](#) for more information

Compliance

EN-ISO 10993-10 (1999-11) / EN 1021-1 and EN 1021-2

Code	Display Name
BTC067500	Closed Seat - 380mm x 425 mm
BTC067900	Closed Seat - 425mm x 425mm

Code	Display Name
BTC067600	Closed Seat - 425mm x 470mm
BTC067300	Closed Seat - 470mm x 470mm