

Pride Litestream Folding Carbon Fibre Power Wheelchair



Effortless, Lightweight and Foldable Travel Companion. With carbon fibre technology, the Litestream[®] Carbon Power Wheelchair is a great travel companion. In just a few quick steps, the Litestream[®] Carbon folds down, making it a super lightweight effortless option to take when travelling in planes, trains and cars. It's a convenient mobility solution with the ability to fold and go.

- Conveniently folds down for transportation.
- Lightweight carbon fibre frame
- Comes with integrated lithium battery backup, [watch our battery maintenance video here.](#)
- Under seat storage
- Turning radius 120 cm
- Joystick controller
- 15 km Range per Charge (varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition.)

Warranty Information

- Frame & Electrical Components: 2 years
- Battery: 1 Year

✓ Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

Compliance

AS 3695.2:2019

Code	Display Name	Total Length	Total Width	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Maximum Speed	Range
PWS654185	Litestream Folding Carbon Fibre Electric Wheelchair	930 mm	577 mm	136 kg	Carbon Fibre	16 kg (exc. Batt.)	430 mm	380 mm	475 mm	475 mm	6 km/h	15 km

