

ROHO - Mosaic



permobil

ROHO

Experience the Value: The re-engineered, cost-effective ROHO MOSAIC cushion is now even better. With a more comfortable fit, you'll enjoy more sitting time while still protecting your skin and soft tissues.

- The innovative design of the MOSAIC cushion allows you to sink comfortably into it, shaping to the unique contours of your body
- MOSAIC provides basic preventative skin and soft tissue protection at an affordable price
- Designed to provide a comfortable, functional and stable sitting environment for patients who are at risk of skin breakdown
- The MOSAIC is also a Dry Floatation cushion but is made from softer plastic than the classic ROHO DF cushions. It is therefore lighter (0.5 kg) but less rugged
- Interconnected air cells disperse pressure and promote stability as the pelvis sinks into them
- The MOSAIC is also excellent as a lounge-chair cushion and can also be used as a shower chair cushion
- All models include: 2 way stretch cover, Hand pump & Repair kit, Operating instructions.

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Compliance

EN-ISO 10993-10 (1999-11) / EN 1021-1 and EN 1021-2

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC586100	ROHO Mosaic - 16 x 16in	2.75 in	143 kg	Low Risk	Air Cells
PTC586200	ROHO Mosaic - 18 x 16in	2.75 in	143 kg	Low Risk	Air Cells
PTC586300	ROHO Mosaic - 18 x 18in	2.75 in	143 kg	Low Risk	Air Cells