

Pride Go-Go Boot Scooter



The Boot Scooter is loaded with design features like one-hand feather-touch disassembly and a convenient drop-in battery box that makes transport and travel worry free. The Boot Scooter delivers these features and more, along with the performance you expect from the first name in travel mobility, making it the ultimate travel scooter value.

- Feather touch disassembly into 5 pieces for easy transportation and storage (refer to image)
- Pride's exclusive black non-scuffing tyres
- Comfortable front-mounted post seat with arm rests for stability
- Easy to remove rechargeable battery
- Soft and easy to adjust tiller
- Flat free tyres
- Front basket standard
- Modular design for easy servicing
- Convenient off-board dual voltage charger that can charge battery pack on-board or off-board

Warranty Information

- 2 Years
- 6 Months Batteries



Contract approved

Support @ Home

Function: Diverse motorised vehicles

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122303122_0105_1_2

Scooter –

Composite:05_418_0105_1_2

Click [here](#) for more information

Compliance

AS3695.2:2013

Code	Display Name	Total Length	Total Width	Safe Working Load	Maximum Speed	Range	Colour
SCT670613	GoGo Boot Scooter - 4 Wheel	1003 mm	495 mm	118 kg	6.4 km/h	11.5 km	Red

