

Pride Go-Go LX with CTS Travel Scooter



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

- Go-Go LX with CTS Suspension (Comfort-Trac Suspension) brings a new level of performance and value to travel mobility
- Comfort Track Suspension (CTS) offers users all the benefits of the Go-Go family, but with the comfort of suspension
- Lightweight mobility scooter
- Disassembles into 5 pieces for easy transportation and storage (refer to image)
- Easy to remove rechargeable battery
- Soft and easy to adjust tiller
- Includes 2 sets of easily changeable colours (Red & Blue)
- Front basket as standard
- Unit Weight: 39.5 kg (without batteries)

Warranty Information

- 2 Years
- 6 Months Batteries



Contract approved

Support @ Home

Function: Diverse motorised vehicles

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122303122_0105_1_2

Scooter –

Composite: 05_418_0105_1_2

Click [here](#) for more information

Compliance

AS3695.2:2013

Code	Display Name	Total Length	Total Width	Safe Working Load	Unit Weight	Colour	Maximum Speed	Range
SCT670612	GoGo LX4 Scooter - 4 Wheel	1040 mm	520 mm	136 kg	39.5 kg (w/out batteries)	Red	6.4 km/h	17 km