

Aspire Vida Folding Manual Wheelchair - Attendant Propelled



- The Aspire Vida is a super lightweight transport wheelchair ideal for quick trips to the shops, appointments and catching up with friends.
- Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support.
- The 135kg safe working load caters to a broad range of users.
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers.
- Contoured and padded seating improves comfort and postural support.
- Attendant operated hand brakes with parking lock.
- Simple and intuitive folding for easy transportation.
- Removable leg rests and folding back for compact storage.
- Rear storage pocket for essentials.
- Available as attendant propelled or self-propelled.
- Anterior Seat Height - 495 mm
- Posterior Seat Height - 445 mm
- Anti-tippers available (WCA848275)
- Lap belt available (WCA848300)

Warranty Information

- 5 Years - Frame
- 1 Year - Fabric & Components



Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3695.1:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Frame Weight	Backrest Height
MWS449850	Aspire Vida Folding Wheelchair - AP - Orange	995 mm	655 mm	930 mm	135 kg	9.6 kg	Orange	420 mm	450 mm	8.8 kg	430 mm
MWS449855	Aspire Vida Folding Wheelchair - AP - Green	995 mm	655 mm	930 mm	135 kg	9.6 kg	Green	420 mm	450 mm	8.8 kg	430 mm