

Aspire Montana HD Lift Recline Chair



- True Bariatric Seating Solution: 300kg and 450kg options available
- Dual Action: Unlike the standard Montana which is a single action, this chair is available in a dual action configuration with independent control of the back, leg rest, and lift movements
- Interchangeable Backrests: Three configurations that can easily be changed due to their velcro fastening with standard (pictured), waterfall, and reinforced lateral support
- Custom Sizing: Capability to offer clients a custom-sized lift chair with the ability to specify the seat width, depth, height and armrest height. This option is available via the custom order form
- Head Cushion: Included as standard (if chosen) for extra comfort and support
- Waterfall and Immersion Backrests are custom and need to be ordered
- Legrest capacity
 - Montana HD (CHP227405): 160 kg
 - Montana HD450 (CHP227410): 140 kg



Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Warranty Information

- 2 Years Mechanism & Electronics
- 2 Years Upholstery
- Limited Lifetime Frame

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Backrest Height
CHP227405	Aspire Montana HD Lift Recline Chair - 254kg	900 mm	1270 mm	254 kg	500 mm	650 mm	810 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Backrest Height
CHP227410	Aspire Montana HD450 Lift Recline Chair - 450kg	1100 mm	1270 mm	450 kg	500 mm	850 mm	810 mm