

## Aspire Shower Chair



- Safe Working Load (SWL) of 200kg tested to AS/NZS39737:2009
- Contoured plastic seat with drainage holes and windows - Provides postural support, improves drainage and reduces skin adhesion
- Height adjustable legs - Easy push button height adjustment for changing needs
- Additional height increments on one leg - Accommodates uneven surfaces
- Raised armrests - Assist the user transferring on and off the chair
- All frames are individually treated for corrosion resistance
- Durable non-slip feet - Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability
- Overall Height: 735 - 910 mm
- Overall Width: 512 mm
- Overall Depth: 570 mm
- Width Between Arms: 470 mm
- Drainage holes on the seat

### Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Washing, bathing and showering

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_093307070\_0103\_1\_2

**Shower Support –**

**Composite:**05\_415\_0103\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTS118005	Shower Chair with Arms - Aluminium	400 mm	440 mm	375 - 550 mm	200 kg	3.1 kg

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTS118000	Shower Chair with Arms - Treated Steel	400 mm	440 mm	375 - 550 mm	200 kg	3.4 kg