

Vogue Lightweight 2 - Seat Walker



This modern seat walker has a clean and crisp frame design, making it a standout when going for a walk or browsing the shops!

LIGHTWEIGHT ALUMINIUM BODY: Improves manoeuvrability and portability.

- Advanced locking breaks
- Unbreakable TPE super soft wheels
- New wider curved backrest
- Minimised storage space with folded width 220 mm
- Locks when folded
- Removable utility bag with 5 kg limit
- Cane holder attachment included
- *Unit weight includes accessories, bag and backstrap
- **Frame weight is minimal functional weight

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Walking, manipulated by both arms

Category: Prescribed

Click [here](#) for more information

Compliance

ISO 11199-2:2005

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Width	Seat To Floor Height	Frame Weight
WAF705310RE	Vogue Lightweight 2 - Red	665 mm	620 mm	785 - 920 mm	150 kg	6.45 kg	Red	465 mm	540 mm	5.3 kg
WAF705310CH	Vogue Lightweight 2 - Champagne	665 mm	620 mm	785 - 920 mm	150 kg	6.45 kg	Champagne	465 mm	540 mm	5.3 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Width	Seat To Floor Height	Frame Weight
WAF705310LE	Vogue Lightweight 2 - Leopard	665 mm	620 mm	785 - 920 mm	150 kg	6.45 kg	Leopard	465 mm	540 mm	5.3 kg