

## Raz CAT Tilt & Recline Shower Commode



- For paediatrics, adolescents or smaller adults, the compact Raz-CAT provides a fit that is better suited to physically smaller clients
- The base and seat frames retain the 18" width of the Raz-AT so that the chair can still be positioned over a toilet
- The depth-adjustable back provides growth capability - essential component for pediatrics



### Support @ Home

**Function:** Washing, bathing and showering

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Seat Depth	Safe Working Load
BTC091800	Raz CAT 460 Tilt & Recline Shower Commode - SA Footrest - Attendant Propelled (as per script)	280 - 400 mm	115 kg