

Aspire ComfiMotion Activ Care Bed



The Activ Care bed combines the best of in-home care functionality, positioning and low height safety with luxurious and contemporary style for enhanced comfort, care and wellbeing.

- **Ultra Low Height:** Minimum platform height of 195 mm reduces falls risk, improving user safety
- **Customisable Safe Entry/Exit Height:** Programmable single button press ensures perfect bed height to maximise independence
- **Under Bed Clearance:** provides clear access for patient lifters and overbed tables
- **Battery Backup:** Provides confidence and maintains independence during power outages
- **Night Assist:** Luminous underbed halo for safe entry and exit
- **Illuminated Handset:** Large button interface is intuitive and highly visible at all times
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief
- **Hi-Lo Functionality:** Assists with manual handling and helps to mitigate fall risks
- **Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV
- **Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea
- **Auto-Regression in the Backrest:** Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution
- **Clinical Positioning:** One touch raises back and legs to alleviate cardiac exertion, improving resting comfort
- **Leg Raise:** Promotes circulation, relieving tired, swollen/ fatigued or painful



Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 60601-2-52:2009

limbs

- **Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button
- **Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back
- **Seated Position:** Assists with bed entry and exit and upright dining
- **Trendelenburg (and Reverse):** Can assist care givers when repositioning a user*
- **Lockable Castors:** Eight heavy duty dual castors
- Headboard and accessories sold separately

*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.

Warranty Information

- 10 Years Limited

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEB046740	Aspire ComfiMotion Activ Care Bed - Portobello Grey - Long Single (Headboard Sep)	2030 mm	920 mm	195 - 700 mm	250 kg
BEB046750	Aspire ComfiMotion Activ Care Bed - Portobello Grey - King Single (Headboard Sep)	2030 mm	1070 mm	195 - 700 mm	250 kg
BEB046760	Aspire ComfiMotion Activ Care Bed - Portobello Grey - Long Double (Headboard Sep)	2030 mm	1370 mm	195 - 700 mm	360 kg
BEB046770	Aspire ComfiMotion Activ Care Bed - Portobello Grey - Queen (Headboard Sep)	2030 mm	1520 mm	195 - 700 mm	360 kg