

Nimbo Lightweight Walker Accessories



Accessories to fit Nimbo Lightweight Walker

- **Seat** - Folds up and out of the way for walking and flips down for sitting. This seat can be added to a Nimbo walker at any time
- **Seat Harness** - Properly positions pelvis and moves with the natural gait pattern of the user, harness is adjustable, available in two sizes
- **Pelvic Stabiliser** - Provides added stability to pelvis for better ambulation. Width adjustable and includes pelvic belt
- **Forearm Platform** - Height, depth and angle adjustable for proper arm support, available in two sizes. (Measurements provided are of the padded plate)
- **Bumper Guard** - Rubber guards protect front of walker and walls from damage
- **Dial-a-speed Tabs** - Allows for variable resistance on rear wheels to control walker speed
- **Outdoor Front Wheels** - Larger 9" non swivel wheels allows single direction movement, sold as a pair, adds 1" to overall height



NDIS

Item Code:

05_180939187_0103_1_2

Walking Aids - Composite:0

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height
PED480900	Nimbo Seat - Tyke (Size 1)	-	-	-
PED480905	Nimbo Seat - Junior (Size 2)	-	-	-
PED480810	Nimbo Seat - Youth (Size 3)	-	-	-
PED480915	Nimbo Seat - Young Adult (Size 4)	-	-	-

Code	Display Name	Total Length	Total Width	Total Height
PED480920	Nimbo Seat - Adult (Size 5)	-	-	-
PED480960	Nimbo Seat Harness Large - Young Adult (Size 4), Adult (Size 5)	-	-	-
PED480965	Nimbo Seat Harness Small - Tyke (Size 1), Junior (Size 2), Youth (Size 3)	-	-	-
PED480950	Nimbo Pelvic Stabilizer	-	-	-
PED480940	Nimbo Forearm Platform Mounting Bracket - Large (Pair)	178 mm	114 mm	89 - 240 mm
PED480945	Nimbo Forearm Platform Mounting Bracket - Small (Pair)	240 mm	127 mm	89 - 240 mm
PED480925	Nimbo Bumper Guard	-	-	-
PED480955	Nimbo Dial-A-Speed Tabs (Pair)	-	-	-
PED480930	Nimbo Outdoors Front Wheels 9" (Pair)	-	-	-