

Aspire Raphael Quattro Lift Recline Chair - Four Action



- Four motor lift recline chair for infinite positioning
- Backlit hand control with integrated USB charging point
- Comes with integrated lithium battery backup, watch our [battery maintenance video here](#).
- Independent adjustable headrest, back, leg & lumbar support functions
- Integrated lift & recline functions for easy sit-to-stand & maximum comfort
- Intuitive headrest support, ideal when watching TV when reclined
- Head & armrest covers for improved longevity
- Part of the Aspire GOLD range of lift recline chairs
- Available in luxurious microsuede and vinyl upholstery options to bring style & class to any environment.
- Legrests extend to provide a longer support under the lower legs, beneficial for taller people

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Compliance

ISO13934-1 / ISO13937-1 / BS6807
/ BS7175

Did you know we can customise your Raphael? - Engage our talented team of Sydney based upholsters to design a unique chair to suit your environment or personality.

Choose from:

- [15 x La Casa is a luxurious artisan leather options](#)
- [15 x Luster Range high performance vinyl options](#)
- [14 x Eastwood Easy Clean Suede options](#)
- [40 x Mazza and Manisa Soft Touch fabric options](#)

Prescribing Therapists: Click here to download our [Clinical Assistance Guide](#)

Code	Display Name	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227290	Aspire Raphael Quattro Powerlift Recline Chair - Small - Bronze	740 mm	1060 mm	130 kg	Bronze Microsuede	480 mm	480 mm	470 mm	470 mm	630 mm
CHP227291	Aspire Raphael Quattro Powerlift Recline Chair - Small - Marina Grey	740 mm	1060 mm	130 kg	Marina Grey Microsuede	480 mm	480 mm	470 mm	470 mm	630 mm
CHP227300	Aspire Raphael Quattro Powerlift Recline Chair - Medium - Bronze	780 mm	1090 mm	130 kg	Bronze Microsuede	500 mm	520 mm	470 mm	470 mm	630 mm
CHP227301	Aspire Raphael Quattro Powerlift Recline Chair - Medium - Marina Grey	780 mm	1090 mm	130 kg	Marina Grey Microsuede	500 mm	520 mm	470 mm	470 mm	680 mm
CHP227302	Aspire Raphael Quattro Powerlift Recline Chair - Medium - Nordic Stone Vinyl	780 mm	1090 mm	130 kg	Nordic Stone Vinyl	500 mm	520 mm	470 mm	470 mm	680 mm
CHP227310	Aspire Raphael Quattro Powerlift Recline Chair - Large - Bronze	810 mm	1080 mm	130 kg	Bronze Microsuede	530 mm	550 mm	500 mm	500 mm	700 mm
CHP227311	Aspire Raphael Quattro Powerlift Recline Chair - Large - Marina Grey	810 mm	1080 mm	130 kg	Marina Grey Microsuede	530 mm	550 mm	500 mm	500 mm	700 mm
CHP227312	Aspire Raphael Quattro Powerlift Recline Chair - Large - Nordic Stone Vinyl	810 mm	1080 mm	130 kg	Nordic Stone Vinyl	530 mm	550 mm	500 mm	500 mm	700 mm
CHP227319	Aspire Raphael Quattro Powerlift Recline Chair - Maxi - Bronze	930 mm	1050 mm	180 kg	Bronze Microsuede	540 mm	650 mm	510 mm	510 mm	650 mm
CHP227322	Aspire Raphael Quattro Powerlift Recline	930 mm	1050 mm	180 kg	Marina Grey	540 mm	650 mm	510 mm	510 mm	650 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
	Chair - Maxi - Marina Grey				Microsuede					