

Aspire Deluxe HD 4 Wheel Scooter - HS898



The HS898 allows for mobility over challenging terrain. Heavy Duty weight capacity of 181 kg, along with a comfortable captain's seat with angle adjustable back support, depth and swivel adjustment on the seat base and angle adjustable tiller. A deluxe heavy duty mobility scooter with hill climbing performance!

- Heavy Duty Hill Climber
- Stylish mag wheels
- LED headlights, brakelights and indicators
- Front and rear wheel suspension
- Rearview mirrors
- Adjustable swing away armrests
- Electro mechanical brakes
- Seat Width 520 mm

✓ Contracted DVA product

✓ Contract approved

Compliance

AS3695.2:2013

Warranty Information

- Scooter Frame: 3 Years
- Motors/Axles/Brakes: 18 Months
- Battery: 6 Months

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Maximum Speed	Range
SCT671290	Aspire Deluxe HD 4 Wheel Scooter - HS898 -	1440 mm	720 mm	1360 mm	181 kg	117 kg (inc.	Black	10 km/h	40 km

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Maximum Speed	Range
	Black					battery)			
SCT671300	Aspire Deluxe HD 4 Wheel Scooter - HS898 - Blue	1440 mm	720 mm	1360 mm	181 kg	117 kg (inc. battery)	Blue	10 km/h	40 km