

Liko HygieneVest (Mod. 55)



- Designed for safe and comfortable transfers to and from the toilet
- Offers almost fully upright sitting, providing functional and safety advantages
- Works well in complicated situations, such as when lifting patients with poor muscle tone
- Vest section firmly supports the upper body during transfer, while a generous opening leaves the lower body free for dressing and undressing
- Available in Polyester

Warranty Information

- 1 Year



NDIS

Item Code:

05_123603166_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011



Code	Display Name	Safe Working Load	Material
LSS394220	HygieneVest Highback Mod 55 Walking Sling - Reinforced Leg Support - Polyester - SMALL	200 kg	Polyester (Loops)
LSS394230	HygieneVest Highback Mod 55 Walking Sling - Reinforced Leg Support - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394240	HygieneVest Highback Mod 55 Walking Sling - Reinforced Leg Support - Polyester - LARGE	200 kg	Polyester (Loops)
LSS394250	HygieneVest Highback Mod 55 Walking Sling - Reinforced Leg Support - Polyester - XLARGE	200 kg	Polyester (Loops)