

## Aspire Deluxe Bedside Commode



- Allows for toileting in the bedroom for users unable to travel between the bedroom and bathroom
- Available in Champagne or Slate Vinyl
- Easy Clean Vinyl Surfaces
- Removable Bowl with Lid & Carry Handle
- Height Adjustable
- Raised Padded Armrests
- Fold Away, Padded Seat
- Non-Slip Tips
- Seat Height: 410 - 565 mm

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Toileting

**Category:** Prescribed

Click [here](#) for more information

### Warranty Information

- 1 Year



### NDIS

**Item Code:**

05\_091203048\_0103\_1\_2

**Toilet Aids – Composite:**

05\_414\_0103\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

| Code      | Display Name                              | Seat Depth | Seat Width | Seat To Floor Height | Safe Working Load | Unit Weight |
|-----------|-------------------------------------------|------------|------------|----------------------|-------------------|-------------|
| BEC045115 | Aspire Deluxe Bedside Commode - Champagne | 450 mm     | 460 mm     | 480 - 630 mm         | 160 kg            | 12 kg       |
| BEC045116 | Aspire Deluxe Bedside Commode - Slate     | 450 mm     | 460 mm     | 480 - 630 mm         | 160 kg            | 12 kg       |

