

Aspire Shower Stool



- Safe Working Load (SWL) of 175 kg - Tested to AS/NZS39737:2009
- Clip on seat attachment - For easy assembly and cleaning
- Height adjustable legs - Easy push button height adjustment for changing needs
- Additional height increments on one leg - Accommodates uneven surfaces
- All frames are individually treated for corrosion resistance
- Raised armrests - Assist the user transferring on and off the chair
- Seat with anterior tilt - Facilitates standing from sitting and provides greater comfort for stiff or painful hips
- Raised armrests - Assist the user transferring on and off the stool
- Durable non-slip feet - Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability
- Available in Aluminium or Steel
- Width Between Arms: 495 mm

Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_093307074_0103_1_2

Shower Support – Composite:

05_415_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTS119005	Shower Stool with Arms - Aluminium	310 mm	375 mm	470 - 655 mm	175 kg	3.1 kg

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTS119000	Shower Stool with Arms - Treated Steel	310 mm	375 mm	470 - 665 mm	175 kg	3.8 kg