

Aspire Midi Deluxe 4 Wheel Scooter - HS520



The ideal Mid-Size four wheel mobility scooter! The HS520 delivers performance with sporty styling. With standard features like a deluxe seat, full suspension and LED Headlights, it will keep you moving.

- Sleek, Sporty Design
- Stylish Mag Wheels
- LED Headlights and Brake lights
- Front and Rear Wheel Suspension
- Rearview Mirror
- Adjustable Swing Away Armrests
- Electro Mechanical Brakes

✓ Contracted DVA product

✓ Contract approved

Compliance

AS3695.2:2013



Warranty Information

- Scooter Frame: 3 Years
- Motors/Axles/Brakes: 18 Months
- Battery: 6 Months

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Range
SCT671210	Aspire Midi Deluxe 4 Wheel Scooter - HS520 - White	1200 mm	640 mm	1155 mm	136 kg	94 kg (inc. battery)	White	39 km