

Ki Mobility Arc Little Wave



- Repositioning - Tilting the occupant has been made easier through an Arc shaped roller guide system which helps translate the Center of Gravity throughout the tilt process
- Attendant Wheel Lock - Easy access to the wheel lock for the attendant
- Drum Brakes
- Folding Recline Back for convenient and fast folding
- Tilt the occupant up to 45°
- Leg extension and chair width have plenty of size option adjustments to accomodate for growing children!

✓ Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
MWC427015	Arc Little Wave Children's Wheelchair (As per script)	300 - 510 mm	250 - 460 mm	370 - 480 mm	75 kg	10 kg