

Aspire XL 4 Wheel Scooter - HS928



This heavy duty scooter is ideal for larger users and for traveling longer distances with a travel range up to 50 km! The sporty style and numerous safety features ensures it is both comfortable and safe, and a real eye catcher. A robust MAXI class scooter that really will go the distance!

- MAXI scooter to cater for larger users
- Stylish mag wheels
- LED headlights and brakelights
- Front and rear wheel suspension
- Rearview mirrors
- Adjustable swing away armrests
- Electro mechanical brakes

Warranty Information

- Scooter Frame: 3 Years
- Motors/Axles/Brakes: 18 Months
- Battery: 6 Months

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Diverse motorised vehicles

Category: Prescribed

Click [here](#) for more information

Compliance

AS3695.2:2013

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Maximum Speed	Range
SCT671310	Aspire XL 4 Wheel Scooter - HS928 - Silver	1600 mm	720 mm	1280 mm	226 kg	178 kg	Silver	10 km/h	50 km