

Aspire Midi Premium 4 Wheel Scooter - HS589



This scooter combines a classic design with smooth, quiet operation as well as the power, travel range and stability needed for outdoor use over varying terrains. A Mid-Size four wheel mobility scooter that's perfect for use outdoors!

- LED headlights and brakelights
- Front and rear wheel suspension
- Rearview mirror
- Adjustable swing away armrests
- Electro mechanical brakes
- Seat Width 455 mm

✓ Contracted DVA product

✓ Contract approved

Compliance

AS3695.2:2013

Warranty Information

- Scooter Frame: 3 Years
- Motors/Axles/Brakes: 18 Months
- Battery: 6 Months



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Maximum Speed	Range
SCT671260	Aspire Midi Premium 4 Wheel Scooter - HS589 - Blue	1290 mm	625 mm	1150 mm	135 kg	100kg (inc. battery)	Blue	10 km/h	40 km
SCT671270	Aspire Midi Premium 4 Wheel Scooter - HS589 - Red	1290 mm	625 mm	1150 mm	135 kg	100kg (inc. battery)	Red	10 km/h	40 km