

Aspire Large Deluxe HD 4 Wheel Scooter - HS828



The HS828 features LED headlights and LCD display with dual side mirrors for that sporty look and performance. A sturdy weight capacity of 160 kg and comfortable and adjustable captain's seat with adjustable tiller. A Large scooter with sporty styling and performance to match.

- Modern design with digital user interface
- Stylish mag wheels
- Headlights, brakelights and indicators
- Front and rear wheel suspension
- LED daytime running headlights
- Rearview mirrors
- Adjustable swing away armrests
- Electro mechanical brakes
- Seat Width 515 mm

Warranty Information

- Scooter Frame: 3 Years
- Motors/Axles/Brakes: 18 Months
- Battery: 6 Months

✓ Contract approved

Support @ Home

Function: Diverse motorised vehicles

Category: Prescribed

Click [here](#) for more information

Compliance

AS3695.2:2013

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Range
SCT671285	Aspire Large Deluxe HD 4 Wheel Scooter - HS828 - Red	1400 mm	680 mm	1325 mm	160kg	124 kg (inc. battery)	Red	40 km