

Aspire Mini Manual-Folding Mobility Scooter - HS268



This compact four wheel mobility scooter is easily transported or stored! Great mobility scooter for smaller, quick trips such as getting around the local shopping mall. As it folds down, the scooter can be stored easily and assembled quickly when needed, making sure you look great on the streets!

- Manual folding via hand lever
- LED headlights and brakelights
- Intuitive user control interface
- Approved for air travel
- Electro Mechanical Brakes
- Fast charging Lithium Ion battery
- Seat Width 419 mm
- Folded Size: 495 (H) x 775 (L) x 510 (W) mm

✓ Contract approved

Compliance

AS3695.2:2013

Warranty Information

- Scooter Frame: 1 Year
- Motors/Axles/Brakes: 1 Year
- Battery: 6 Months

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Maximum Speed	Range
SCT671135	Aspire Mini Manual-Folding Mobility Scooter -	950 mm	500 mm	960 mm	115 kg	30.6 kg (inc.	White	8 km/h	20km

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Maximum Speed	Range
	HS268 - White					battery)			