

Ceiling Hoist - Liko - Likorall 200



- The Likorall 200 is an economical fixed ceiling hoist needed to manage everyday lifting situations
- The patented safety drum with speed limiter facilitates electrical emergency lowering
- Ideal for floor lifts in most environments
- Can be equipped with a special carriage for planned transfers between different rooms
- Ergonomically designed hand control
- Quick sling bar change thanks to Q-link and Quick-release Hook
- Simple wall charging via hand control

✓ Contracted DVA product

✓ Contract approved

Compliance

AS/NZS ISO 10535:2011

Warranty Information



- Frame: 10 Years
- Electrical Parts: 2 Years

Code	Display Name	Total Height	Safe Working Load
LIC359490	Likorall 200 - Fixed Ceiling Hoist	Up to 1800 mm	200 kg