

WHILL MODEL C2 / F - ACCESSORIES



Accessories for your Whill C2 and Whill F power wheelchair

- Large under seat Storage basket securely stores your belongings and groceries (Model C2 only).
- Side storage bag ideal for personal belongings such as phone and wallet (Model C2 only).
- Dust Cover: Protects the device from dust and dirt when stored outdoors (Do not use for rain protection, cover not waterproof)(Model C2 only).
- Lap Belt: Firmly stabilises the sitting position even on slopes.
- Smart Key: It provides you with the ability to lock and unlock your Whill from a distance, similar to a car's key fob.
- Battery: The 2.7kg lightweight battery can also be used as a spare battery.
- Battery Charger: A great compliment to a spare battery, it helps let you go from one fully charged battery to the next fully charged battery to keep you moving.
- Whill C2 Power Wheelchair standard inclusions: Lap Belt, Smart Key, Battery & Battery Charger.
- Whill F Power Wheelchair standard inclusions: Battery & Battery Charger.



Contract approved

Support @ Home

Function: Wheelchair accessories

Category: Prescribed

Click [here](#) for more information

Code	Display Name
WCA870260	Whill Model C2 - Dust Cover
WCA870400	Whill Model C2 / F - Lithium Ion Battery
WCA870410	Whill Model C2 / F - Charger

Code	Display Name
WCA870420	Whill Model C2 / F - Lap Belt
WCA870430	Whill Model C2 - Side Storage Bag
WCA870440	Whill Model C2 - Storage Basket
WCA870450	Whill Model C2 - Walking Stick Holder
WCA870710	Whill Model F Accessory - Walking Stick Holder
WCA870700	Whill Model F Accessory - Smart Key