

Frontier V6 - Compact 40



magic
MobilityTM

The baby of the family, the Compact 40 is ideal for negotiating tight turns. Small but perfectly formed, the only thing that we've reduced is the size. Still the same frame, motors and seating. We just made it smaller.

✓ Contract approved

- The Magic Mobility Frontier V6 all-terrain mid-wheel drive power wheelchair provides tremendous off-road and curb-climbing capabilities
- Its low pressure tires provide a smoother ride and give the Frontier V6 increased climbing ability
- POWER SEAT ELEVATE for children: Being at eye level to their fellow classmates and friends through the use of the power lift function can help reduce the negative stigma that may affect young wheelchair users. This can result in more enriching positive experiences for children.
- Greatest off-road capability in the Frontier Range
- Great for curb climbing
- Exceptionally smooth ride with low pressure tyres
- Crash test compliant tie down points
- TAKE CONTROL - A full range of advanced control systems and environmental controls are available. Easy to program and easy to use
- POWER OPTIONS - Consider a seat lift to help you reach, transfer and face life head on, or power tilt for improved manoeuvring over obstacles.
- UP CLOSE AND PERSONAL - Side transfers aren't the only reason you need to flip your armrest out of the way
- Gel Batteries
- Drive mode: Mid wheel drive
- Motors: 68Nm Torque, 800W, Gear In Line
- Environment: City streets and tight indoors
- Full lighting packages
- Rehab positioning:

- Power options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swing away, power or gas and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available

Code	Display Name	Total Length	Total Width	Safe Working Load	Seat To Floor Height	Seat To Footrest	Maximum Speed
PWC655050	Frontier V6 - Compact 40 (as per script)	910 mm	655 mm	182 kg	435 mm	435 mm	10 km/h