

JAY - Duo



JAY[®]

- Designed for the client who has symmetrical posture needing mild to moderate physical support
- The lightweight aluminium pre-contoured solid base attaches to chair rails effectively replacing the sling seat upholstery
- The aluminium base in combination with a soft foam & fluid cushion provides an anatomically correct, stable base of support that promotes postural stability
- The pre-contoured foam base of the Jay Duo provides mild lateral and forward/rearward stability and provides a base of support that promotes foot propelling or independent transfers
- The 'Jay Flow' fluid pad reduces the effects of shearing on the ischial tuberosities
- Most suitable for the client at moderate risk of skin breakdown who can ideally perform independent weight-shift
- All Models Includes 2 way stretch cover

✓ Contracted DVA product

✓ Contract approved

NDIS

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Warranty Information

- 2 Years Cushion
- 6 Months Cover

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC569200	JAY Duo - 16 x 16in	2 in	113 kg	Medium Risk	Foam/Fluid Pad
PTC569300	JAY Duo - 16 x 18in	2 in	113 kg	Medium Risk	Foam/Fluid Pad

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC569500	JAY Duo - 17 x 18in	2 in	113 kg	Medium Risk	Foam/Fluid Pad
PTC569000	JAY Duo - 18 x 16in	2 in	113 kg	Medium Risk	Foam/Fluid Pad
PTC569100	JAY Duo - 18 x 18in	2 in	113 kg	Medium Risk	Foam/Fluid Pad
PTC569400	JAY Duo - 20 x 18in	2 in	113 kg	Medium Risk	Foam/Fluid Pad