

Vicair Active O2 Cushion



- The Active O2 offers the unique combination of ergonomically shaped foam and Vicair SmartCells
- The four rear SmartCell compartments provide a high level of comfort and skin protection, and reduce pressure on the sitting bones and tailbone area
- The foam front makes transfers easy and feels steady and improves wheelchair handling
- The foam inlay is protected by an easy-to clean fluid-resistant liner. The rest of the cushion is 100% machine washable and breathable for optimal hygiene and perfect temperature & moisture regulation
- Perfect for active wheelchair users, who transfer frequently
- The Active O2 has an easy accessible storage pouch for precious belongings (keys, phone, wallet)
- Trochanter support firmer density to provide lateral stability - redistributes pressure to help create an ischial well
- SmartCell - Foam barrier limits sliding forward
- Ergonomically-shaped foam for easy transfers
- Ischial well - Two compartments for superior pressure reduction in the IT-coccyx area & extra positioning possibilities



Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTC609460	Vicair Active O2 Cushion 6 - 18x16in	16 in	18 in	6 cm	200 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTC609470	Vicair Active O2 Cushion 6 - 18x18in	18 in	18 in	6 cm	200 kg
PTC609480	Vicair Active O2 Cushion 6 - 18x20in	20 in	18 in	6 cm	200 kg
PTC609490	Vicair Active O2 Cushion 6 - 20x16in	16 in	20 in	6 cm	200 kg
PTC609500	Vicair Active O2 Cushion 6 - 20x18in	18 in	20 in	6 cm	200 kg
PTC609510	Vicair Active O2 Cushion 6 - 20x20in	20 in	20 in	6 cm	200 kg