

Pride Pathrider 10 Scooter



- The sleek, sporty Pathrider 10 scooter delivers high performance operation
- Easy drive tiller with wrap-around handles
- Pride's exclusive stylish, light weight, low profile, non-scuffing, puncture-proof wheels
- New back-lit LED battery meter for increased accuracy and visibility in any setting
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight, folded plastic seat with viscoelastic foam inserts for increased comfort
- Includes pinchless hinges & eight-position seat rotation
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port for standard dual voltage off-board charger
- Front basket standard
- Unit Weight: 63.5 kg (without batteries)

✓ Contracted DVA product

✓ Contract approved

NDIS

Item Code:

05_122303122_0105_1_2

Scooter –

Composite:05_418_0105_1_2

Click [here](#) for more information

Compliance

AS3695.2:2013

Warranty Information

- 2 Years
- 1 Year Batteries

Code	Display Name	Total Length	Total Width	Safe Working Load	Unit Weight	Maximum Speed	Range
SCT670615	Pathrider 10 Mobility Scooter - 4 Wheel (inc Rear bag + Cane	1190 mm	565 mm	159 kg	63.5 kg (w/out	8.5 km/h	25 km

Code	Display Name	Total Length	Total Width	Safe Working Load	Unit Weight	Maximum Speed	Range
	Holder)					batteries)	