

Forearm Crutches



- Lightweight and durable forearm crutches to suit a large range of users
- Cuff provides a secure fitting to the forearm
- Adjustable height between cuff and handle & handle and ground
- Ferrule: 22 mm
- Product Weight: 600 g (each)
- User Height: WAC692000 (Youth Small) = 1230-1640mm, WAC691400 (Adult Medium) = 1480-1860mm, WAC691800 (Adult Tall) = 1740-2040mm

Warranty Information

- 1 Year



✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Walking, manipulated by one arm

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_120309088_0105_1_2

Walking Aids -

Composite: 05_416_0105_1_2

Click [here](#) for more information

Compliance

ISO 11199-2:2005

Code	Display Name	Total Height	Safe Working Load
WAC692000	Forearm Crutches Youth/Small	530 - 785 mm	180 kg

Code	Display Name	Total Height	Safe Working Load
WAC691400	Forearm Crutches Adult/Medium	635 - 910 mm	180 kg
WAC691800	Forearm Crutches Adult Tall/Large	820 - 1075 mm	180 kg