

Liko Viking M Lifter



- Ideal for meeting the lifting needs in small spaces
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales
- Lifting Height 420 - 1790 mm

Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years



Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393190	Liko Viking Lifter - MEDIUM	1230 mm	690 - 1110 mm	1340 - 2020 mm	205 kg	Aluminium
LSS393200	Liko Viking - M - Leg Protector Set -	-	-	-	-	-