

Pride Pathrider 10 DX Scooter



- High performance sleek and sporty mobility scooter
- Low profile, puncture proof tyres
- Comfortable front - mounted post seat with arm rests for stability
- Rugged front and rear suspension
- Long lasting high intensity LED headlight, tail light and indicators
- Back-lit battery gauge on monitor
- Wrap around easy drive tiller
- Tiller mounted battery charger port for off board charging
- Front basket as standard
- Unit weight: 63.5 kg (without batteries)

Warranty Information

- 2 Years
- 1 Year Batteries

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Diverse motorised vehicles

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122303122_0105_1_2

Scooter –

Composite:05_418_0105_1_2

Click [here](#) for more information

Compliance

AS3695.2:2013

Code	Display Name	Total Length	Total Width	Safe Working Load	Unit Weight	Maximum Speed	Range
SCT670620	Pathrider 10 Deluxe Mobility Scooter - 4 Wheel (inc Rear bag +	1193 mm	562 mm	158 kg	85 kg (w/out	9 km/h	25 km

Code	Display Name	Total Length	Total Width	Safe Working Load	Unit Weight	Maximum Speed	Range
	Cane Holder)					batteries)	