

Aspire Overbed Table - Laminate Flat Top



- Lightweight mobile table positions overbed or chair
- Can be used for eating, drawing or other activities
- Flat top ideal for hospital, facility or home use
- Fully height adjustable with a fixing knob for additional security
- C base shape fitted with four castors enables close positioning to the bedside or chair
- Powder coated mild steel frame construction
- Cream standard laminate top with vinyl wrap around the edge, or thermoform cream top which is heat resistant
- Available in White & Beech Laminate coloured tops

Warranty Information

- 1 Year



✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Tables

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_180315322_0103_1_2

Minor Bed – Composite:

05_428_0103_1_2

Click [here](#) for more information

Compliance

ISO7172

Code	Display Name	Total Length	Total Width	Total Height
BEA006700	Aspire Overbed Table - Laminate Flat Top - White	430 mm	800 mm	730 - 1095 mm
BEA006701	Aspire Overbed Table - Laminate Flat Top - Beech	430 mm	800 mm	730 - 1095 mm

