

Whill Model C2 Power Wheelchair



The WHILL Model C2 is a revolutionary approach to power mobility that encompasses all the exciting features of the C. It uniquely boosts user confidence and evokes an active, can-do lifestyle. Its premium performance and superb engineering provide versatility and reliability to experience all-day adventures.

- Mobile phone app compatibility: access the Whill phone app and operate your chair from anywhere, on all Apple & Android devices.
- Improved user interface: ensures all functional controls are in one place, allowing for simple chair control by the user.
- Rear tail lights: improves safety for users in the chair when crossing roads or in high density areas.
- Rear wheel suspension: Absorbs the impact when driving over bump roads or obstacles.
- Up to a 45% tighter turning radius than tiller operated mobility devices.
- Patented front omni-wheels for tight turning and superb terrain coverage.
- Equipped with two powerful motors and engineered for superb stability.
- Electromagnetic brakes provide safe stopping, even on inclines.
- An anti-sway system keeps you driving straight on side slopes.
- Quickly disassembles into three lightweight components for easy transport or storage.
- Charge your mobile devices with a built-in USB port.
- **Standard inclusions:** lap belt, smart key, black side covers, battery and battery charger.

Warranty Information

- Base & Seat Frame: 5 Years
- Electrical & Mechanical Components: 2 years



Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

Compliance

AS 3695.2:2019

- Batteries & Wheels: 6 months

Code	Display Name	Total Length	Seat Width	Total Width	Safe Working Load	Maximum Speed	Range
PWS655400	Whill Model C2 Power Wheelchair - 450mm	985 mm	450 mm	604 mm	136kg	8 km/h	18 km
PWS655410	Whill Model C2 Power Wheelchair - 500mm	985 mm	500 mm	654 mm	136kg	8 km/h	18 km