

Pride Pathrider 140 XL Scooter



- Top of the range sleek and sporty powered mobility scooter
- Large 13" pneumatic tyres to handle rugged outdoor terrain
- Luxurious suspension seat with arm rests, head rest and lumbar support for stability
- Front and rear suspension and sealed hydraulic brakes
- Automatic full lighting system including directional adjustment. Long lasting high intensity LED kerb light, headlight, tail light and indicators
- Back lit battery gauge
- Wrap around easy drive tiller
- Tiller mounted battery charger port for off board charging
- Unit Weight: 118 kg

Warranty Information

- 2 Years
- 1 Year Batteries

✓ Contracted DVA product

✓ Contract approved

NDIS

Item Code:

05_122303122_0105_1_2

Scooter –

Composite:05_418_0105_1_2

Click [here](#) for more information

Compliance

AS3695.2:2013

Code	Display Name	Total Length	Total Width	Safe Working Load	Maximum Speed	Range	Colour
SCT670578	Pathrider 140 XL Scooter - Yellow - 4 Wheel	1428 mm	692 mm	181 kg	10 km/h	35 km	Yellow
SCT670576	Pathrider 140 XL Scooter - Red - 4 Wheel	1428 mm	692 mm	181 kg	10 km/h	35 km	Red
SCT670577	Pathrider 140 XL Scooter - Silver - 4 Wheel	1428 mm	692 mm	181 kg	10 km/h	35 km	Silver