

Aspire Low Back Classic Day Chair



- A versatile, low back orthopaedic chair
- Contoured Backrest – Provides postural support
- Softly padded armrests, seat and back upholstery are comfortable and supportive
- Non-Slip Rubber Tips – Provide stability whilst sitting or transferring
- Easy clean vinyl surfaces
- Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs
- Lightweight Design – lightweight yet sturdy design for easy manoeuvrability
- Powder Coated Steel Frame – Easy to clean, strong and durable
- Available in either Slate or Champagne Vinyl

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

NDIS

Item Code:

05_180915811_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 4688.2:2000; AS/NZS

4688.3:2000

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Colour
CHP208005	Aspire Low Back Classic Day Chair - Slate Vinyl	440 mm	460 mm	440 mm	430 - 605 mm	160 kg	Slate Vinyl