

Aspire Matisse Quattro Lift Recline Chair



- **Independent Adjustable Headrest, Back, Leg and Lumbar** support functions
- **Gentle Curve in Seating Surface** for an improved seating experience
- **Large Side Storage Pockets**, ideal for books and newspapers
- **Integrated Footrest Extension** suits most user heights
- **Integrated Lift and Recline** functions for easy sit-to-stand
- **Backlit Hand Control** with integrated USB charging point
- **Adjustable Headrest**, ideal for watching TV whilst reclined
- **Lithium Ion Battery Backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day
- **4 x Swivel Castors** allow the chair to be moved around with ease

✓ Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Warranty Information

- Frame: 5 Years
- Components*: 2 Years

* Components: Handset, upholstery and electronics.



Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227560	Matisse Quattro Lift Recline Chair - Large	965 mm	1065 mm	159 kg	545 mm	570 mm	500 mm	500 mm	730 mm