


Zippie Iris



- The Zippie IRIS® rotation-in-space wheelchair offers a superior combination of positioning, maneuverability, and increased portability with the XLOCK® folding crossbrace option.
- Cutting-edge rotation-in-space technology designed just for kids
- This tilting wheelchair uses patented rotation-in-space technology to rotate the seat frame around the user's center of gravity.
- This rotational system provides the shortest possible wheelbase without compromising stability to ensure a smooth and controlled rotation.
- The Zippie IRIS also comes with the True Fit Program, a complete growth package for free.

 Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122203115_0105_1_2

Wheelchair Manual –

Composite:05_417_0105_1_2

Click [here](#) for more information

Compliance

AS 3696-1:2008/AS

3696-3:2008/AS 3696-5:1989/AS

3696-8:1998/AS 3695:1992

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
MWC427500	Zippie Iris	36 - 56 cm	25 - 48 cm	32 - 50 cm	125 kg	17 kg