

## Extreme X8 - 4x4



**magic**  
Mobility

- The Extreme X8 has a top speed of 10 km/h (6.2 mph)
- Articulating front frame with all terrain low pressure knobby tyres, which provides an exceptionally smooth and jolt-free ride
- Comes in five standard colours and three specialised prints
- Four independent motors to make kerbs and small steps easy to climb
- Splash resistant (not waterproof)

✓ Contract approved

### Support @ Home

**Function:** Powered wheelchairs

**Category:** Prescribed

Click [here](#) for more information

| Code      | Display Name                       | Total Length | Total Width | Safe Working Load | Seat To Floor Height | Seat To Footrest | Maximum Speed |
|-----------|------------------------------------|--------------|-------------|-------------------|----------------------|------------------|---------------|
| PWC655000 | Extreme X8 - 4x4 (Base Price only) | 1030 mm      | 710 mm      | 182 kg            | 470 mm               | 470 mm           | 10 kph        |