

Liko Original Highback Sling



- Can be used to lift most patients
- Good choice for lifting from a horizontal position
- Provides support for the entire back and head
- Enables a semi-reclined sitting posture

Warranty Information

- 1 Year



Contract approved

NDIS

Item Code:

05_123603164_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS367854	Original High Back MOD 200/210 Net Poly - MEDIUM	200 kg	Net Polyester
LSS367856	Original High Back MOD 200/210 - Net Poly - XLARGE	500 kg	Net Polyester