

Malvestio Brera RX Patient Trolley



- Foot operated hydraulic lift, tilt, back raise, Trendelenburg, and reverse Trendelenburg functions
- Full length sliding X-ray deck ensures complete spinal coverage
- Comfortable 100mm foam mattress with multi-stretch cover (purchased separately)
- Quick Release CPR with access from both sides for backrest and knee break ensures the trolley is emergency ready at all times
- Collapsible vertical side rails fold away below the patient surface
- Fold-away carer push handles to maximise workspace
- 200mm castors minimise rolling resistance with central and directional locking functions
- Molded base cover facilitates cleaning and supports infection control processes
- Large corner bumpers minimise wear and tear on equipment and facilities
- Frame is finished with Biocote, a technologically advanced surface coating that promotes infection control via its anti-bacterial properties
- Backrest angle 0 to 90 degrees
- Trendelenburg -17 to 17 degrees
- **Optional Accessories**
 - IV Pole (TSE674130)
 - Oxygen Bottle Holder (TSE674110)
 - Folding Monitor Tray (TSE674120)



Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
TSE674300	Malvestio Brera RX Patient Trolley	2100 mm	830 mm	560 - 880 mm	270 kg