

Progeo - Easy Tilt



permobil

Progeo

- From Progeo's decade of experience with tilting wheelchairs we have designed Easy Tilt.
- We have set ourselves important objectives such as an excellent quality / price ratio, greater strength and rigidity of the frame, wide adjustments, ease to control tilting movement and finally a reduced size.



Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Backrest Height
MWC451230	Progeo - Easy Tilt (as per script)	900 mm	125 kg	17.5 kg	350 - 500 mm	330 - 480 mm	420 - 540 mm