

Liko Viking L Lifter



- An everyday flexible lifting solution
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales
- Lifting Height 460 - 1790 mm

Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years



Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393180	Liko Viking Lifter - LARGE	1360 mm	695 - 1070 mm	1390 - 2030 mm	250 kg	Aluminium
LSS393210	Liko Viking L - Leg Protector Set -	-	-	-	-	-