

SitnStand - Classic



- Give your life a lift everywhere
- Designed to fit any chair comfortably
- Easy to use, safe and stable
- Weight allowance up to 181kg
- Washable, easy to remove seat cover
- Lift weight 4kg
- Designed for prolonged sitting/easy standing
- Rechargeable battery, no cords.
- Battery life 4-7 days

Warranty Information

- 2 years

✓ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Safe Working Load
LSS403710	SitnStand Lift Assist - Classic	121 kg