

SmartDrive MX2 + PushTracker E2



- **Be Free With SmartDrive MX2+:** The SmartDrive gives you freedom to do more. It moves with you and it's so lightweight, you don't even know it's there. You can spin in place, be in a wheelie, or hop off curbs. It is compatible with your active lifestyle.
 - Ergonomic design quickly attaches to a wide variety of manual wheelchairs
 - Lightest rear power assist in the world
 - OmniWheel design allows SmartDrive to start from any position maintaining optimal traction and maneuverability
- **More connectivity and durability, with a set-up that takes just minutes with PushTracker E2:** The PushTracker E2 makes using your SmartDrive easier than ever and was developed with connectivity and durability in mind. It's been rigorously tested and proven to be consistently reliable and durable. No Wi-Fi or phone required to get started with the SmartDrive MX2+ and PushTracker apps pre-installed. With gesture-activated controls you can start your SmartDrive with a simple tap of your hand—no pushing or reaching necessary.
 - Touch screen technology
 - 8x more reliable Bluetooth® connection
 - Quickly change all performance settings

Warranty Information

- 2 Years



Contract approved

Support @ Home

Function: Wheelchair accessories

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122212120_0105_1_2

Wheelchair Manual – Composite:

Wheelchair Powered –

Composite:05_419_0105_1_2

Click [here](#) for more information

Compliance

EN-ISO 10993-10 (1999-11) / EN
1021-1 and EN 1021-2

| Code | Display Name | Total Length | Total Width | Total Height | Unit Weight |
|-----------|----------------|--------------|-------------|--------------|-------------|
| WCA858508 | SmartDrive MX2 | 389 mm | 141 mm | 242 mm | 5.7 kg |
| WCA858513 | PushTracker E2 | - | - | - | - |