

## PAWS Power Assisted Wheelchair Systems- Tourer



- PAWS Tourer has the longest wheelbase so while the big wheel overcomes roots and uneven terrain, it also makes a smooth ride while touring.
- PAWS Tourer is the true off-road stager. With a 20" wheel and 4" fat tyre, you will easily get to places that would be difficult for other models.
- Designed primarily for off-road purposes, it is not afraid of sand, gravel, stones. Get your PAWS in the dirt!
- Optional Tetra Controls allow you to drive easily, even with diminished upper extremity strength and control.
- With PAWS your manual wheelchair becomes an exciting vehicle, allowing you to be independent when getting to work, going shopping, or on a trip.
- It combines ergonomics and functionality with good looks.
- Every PAWS creates its own story of freedom – meet family and friends, participate in more activities together or go for a solo adventure.
- When choosing PAWS you can customise the set up with different configurations of clamps, steering, wheels, and accessories for optimal use and fit – especially beneficial for quadriplec/tetraplegic users.
- With Auto Lift, you don't need to use upper body strength to connect PAWS with the wheelchair. Auto Lift will automatically connect you to driving position.
- PAWS comes with a variety of speed modes (\*maximum speed is 10 km/h).
- PAWS has an IP54 waterproof rating so is not afraid of rain, puddles and moisture. Go on, get PAWS out!
- The Lithium Ion battery only needs up to 5 hours to be fully charged.
- On a full charge you can reach a range of up to 40 km (depending on terrain, driver weight and model chosen).
- An IATA compliant "flying" battery is also available for when you take your PAWS on a flight.
- With Automatic Easy Clamps just roll your wheelchair into the clamps and



Contract approved

### Support @ Home

**Function:** Wheelchair accessories

**Category:** Prescribed

Click [here](#) for more information

press the button for automatic clamping and lifting to the driving position.  
From Start to go in 12 seconds!

Code	Display Name
WCA872150	PAWS Tourer - 20inch - Auto Mounting - Quad/Tetra Hand Control
WCA872140	PAWS Tourer - 20inch - Auto Mounting - Standard Hand Control